

# MAD Microfocus Macromolecular Crystallography at Diamond Light Source

Robin L. Owen, Danny Axford, Jun Aishma,  
Emma Shepherd, Gwyndaf Evans  
Diamond Light Source, Didcot,  
OX11 0DE, United Kingdom


## Introduction

Increasingly, a critical restraint on the application of X-ray diffraction has been the difficulty in obtaining crystals of sufficient size (50-100 microns) and order to generate the high quality diffraction patterns required for the detailed analysis. This is notably true for membrane proteins and multi-protein complexes, which play a central role in cellular operation. A consequence of the drive to increase the rate of crystal automation of protein expression, screening and crystallisation is the growing dependence on small amounts of material for the structural determination. A solution to this problem is to reduce the beam size and preferentially illuminate smaller crystals or more perfect regions.

## Key Features

Diamond : 3GeV, 300 mA, 2.7 nm.rad  
Source : 2m long 21mm period undulator (> 5mm gap)  
Canted ID angle : 1.5 mrad outboard  
Energy range : 6.5-25 keV  
Energy resolution :  $\sim 2 \times 10^{-4}$  ( $\Delta E/E$  of Si 111)  
Flux @ 12 keV at sample :  $> 10^{12}$  ph/s into focal spot  
Beam size at sample (focused) :  $< 5\mu\text{m}$  (h)  $\times$   $5\mu\text{m}$  (v), FWHM  
Beam diverg. at sample (h  $\times$  v) :  $< 2.0\text{mrad} \times 0.3\text{mrad}$   
Harmonic rejection :  $10^{-4}$


## Beamline design



## User operation

### Experimental Control


Users use the Generic Data Acquisition (GDA) graphical interface to run their experiments in their entirety. This includes upload of sample data from ISPyB, sample centering, interface to robot, and initiation of experiment. Through GDA users have full control of the beamline, with changes in beamsize at the sample and energy changes easily accessible. In addition tools are being developed to allow users to exploit the opportunities made available by a microfocus beamline.



### Grid Scan

This tool allows a search grid to be stretched and dragged over loop/region of interest.


- Beam can be focused/defocused for fine/coarse scan
- Search for small crystals or more perfect regions of larger crystals



### Line Scan

Start and end points for data collection are defined on the crystal via the gui allowing helical data collection.

- Unexposed crystal is continually introduced into the beam during data collection



The second stage of I24 is an ultra-stable microfocusing stage consisting of K-B mirrors, beam conditioning and diagnostics unit, sample stage and goniometer mounted on a common support structure. No feedback will be used within this shorter stage since residual instabilities in the secondary source will be further demagnified to a negligible level.

During the commissioning phase, users have access to two optical arrangements, 'mini' and 'maxi'. The mini configuration results in a beam of  $\sim 9\mu\text{m}$  (h)  $\times$   $9\mu\text{m}$  (v) FWHM at the sample (below). The maxi beam is roughly  $40\mu\text{m}$  (h)  $\times$   $50\mu\text{m}$  (v) FWHM at the sample and is mostly commonly used during the grid scan (details to right).

