

<b>Subject:</b>	<b>SDL End Station Laser Area Interlock Test</b>		
<b>Number:</b>	LS-PPS-0040	<b>Revision:</b>	A
		<b>Effective:</b>	4/13/04
			<b>Page 1 of 6</b>

Prepared/ Approved By:	M. Buckley	Prepared/ Approved By:	S. Buda	Approved By:	
------------------------	------------	------------------------	---------	--------------	--

\*Approval signatures on file with master copy.

[Revision/Periodic Review Log](#)

Test Reason:	Test Result:	<input type="checkbox"/> Passed	<input type="checkbox"/> Failed
	Test Type:	<input type="checkbox"/> Full	<input type="checkbox"/> Partial
Test Date:	Start Time:	Finish Time:	
Tester 1:	Assistant 1:		
Tester 2:	Assistant 2:		

#### Preparation/Prerequisites:

- ☐ Inform the experimenter that a test of the laser interlock is being performed.
- ☐ A qualified SDL laser operator must be available to assist in placing Safeguards (Beam-Dumps) at the output of each laser system as indicated below and to place the laser systems in a "Ready-Mode" where the power supplies will be on with no laser output.

Place beam-dump/safeguards in place and Red-Tag each laser output:

- ☐ Ti-Sapphire laser (Titanium-Sapphire) - located in the Main Laser Hutch (MLH)
- ☐ Nd: YAG laser (neodymium:yttrium-aluminum-garnet) - Main Laser Hutch (MLH)
- ☐ Excimer laser - End Station Laser Area (ESLA)
- ☐ Nd: YAG laser (neodymium:yttrium-aluminum-garnet) - End Station Laser Area (ESLA)

#### 1. **[Search Sequence]**

Press the Start Search button located inside the ESLA.

The start search indicator comes ON

Exit the ESLA, close the door, and press the Complete Search Button.

Observe the steady audible warning is a minimum of 13 - 15 seconds.

Both "Laser Interlocked" signs comes ON

#### 2. Open the door

Intermittent audible alarm sounds for 5 seconds.

Both "Laser Interlocked" signs turn OFF

Interlock fault indicator comes ON

Attempt to secure the ESLA

Cannot complete search - Laser interlock does not activate

Reset Interlock fault

#### 3. Press the start search button. Exit the ESLA, close the door and wait.

The Start Search indicator turns OFF in 40 seconds or less.

Press Complete Search Button.

Interlock does not activate

<b>Subject:</b>	<b>SDL End Station Laser Area Interlock Test</b>		
<b>Number:</b>	LS-PPS-0040	<b>Revision:</b>	A
		<b>Effective:</b>	4/13/04
			<b>Page 2 of 6</b>

4. Press the Start Search button. Exit the ESLA and close the door. Press the Complete Search Button. Open the door and close the door within the time-out period and press the Complete Search Button.

The interlock does not activate

\_\_\_\_\_

5. Search and interlock the ESLA. Press the interlock OFF button

The laser interlock drops out

\_\_\_\_\_

Open the door.

No audible warning sounds.

\_\_\_\_\_

Interlock fault indicator does not come ON

\_\_\_\_\_

6. At this time the ESLA should not be secured.

Verify the Faults are cleared on the two YAG power supplies in the Main Laser Hutch (MLH)

\_\_\_\_\_

Verify Faults are cleared on the Ti-Sapphire power supply laser power supply in the Main Laser Hutch (MLH)

\_\_\_\_\_

Verify Faults are cleared on the Excimer power supply in the End-Station Laser Area (ESLA)

\_\_\_\_\_

Verify Faults are cleared on the (ESLH)YAG power supply in the End-Station Laser Area (ESLA)

\_\_\_\_\_

The ESLA laser shutter is closed.

\_\_\_\_\_

7. Attempt to open the ESLA shutter.

The ESLA shutter does not Open.

\_\_\_\_\_

8. 8.1) Secure the ESLA and MLH for laser operation (Assistant remains inside).

8.2) Request the laser operator to turn ON all of the laser power supplies that will be involved with this test and achieve a 'laser ready' state where no laser beam will be transmitted.

8.3) Open the ESLA laser shutter.

The ESLA shutter opens

\_\_\_\_\_

A short audible sound followed by a 2-second audible signal is heard

\_\_\_\_\_

The 3 power supplies are ON in the MLH

\_\_\_\_\_

The Excimer power supply in ESLA is ON

\_\_\_\_\_

The YAG power supply in ESLA is ON

\_\_\_\_\_

<b>Subject:</b>	<b>SDL End Station Laser Area Interlock Test</b>		
<b>Number:</b>	LS-PPS-0040	<b>Revision:</b>	A
		<b>Effective:</b>	4/13/04
			<b>Page 3 of 6</b>

9. **[Pass Through]**

Press the Internal PASS-Through button. Open the door and exit the ESLA.

The PASS Through button indicator comes

ON

The shutter remains open

The 3 power supplies remain ON in the MLH

The Excimer power supply in ESLA remains

ON

The YAG power supply in ESLA remains ON

Close the door within the timeout period.

9 cont. Press the External PASS-Through button.

The PASS Through button indicator goes out

The shutter remains open

The 3 power supplies remain ON in the MLH

The Excimer power supply in ESLA remains

ON

The YAG power supply in ESLA remains

ON

Laser security remains active in ESLA

10. Press the External PASS-Through button. Open the door and enter the ESLA.

The PASS Through button indicator comes

ON

The shutter remains open

The 3 power supplies remain ON in the MLH

The Excimer power supply in ESLA remains

ON

The YAG power supply in ESLA remains ON

Close the door within the timeout period.

Press the Internal PASS-Through button.

The PASS Through button indicator goes out

The shutter remains open

The 3 power supplies remain ON in the MLH

The Excimer power supply in ESLA remains

ON

The YAG power supply in ESLA remains ON

Laser security remains active

11. **[ESLA Laser Shutter Operation]** Cycle the shutter Open/Close switch.

The shutter opens and then closes with out any hesitation.

<b>Subject:</b>	<b>SDL End Station Laser Area Interlock Test</b>		
<b>Number:</b>	LS-PPS-0040	<b>Revision:</b>	A
		<b>Effective:</b>	4/13/04
			<b>Page 4 of 6</b>

12. Open the laser shutter. Press the External PASS-Through button,  
Open the door, and begin timing.

After no more then 45 seconds the following occurs:

An audible warning sounds for 5 seconds

The interlock fault LED comes ON

Laser Interlock drops out

The shutter closes

The 3 power supplies remain ON in the MLH

The Excimer power supply in ESLA Turns  
OFF

The YAG power supply in ESLA Turns OFF

Interlock fault will only reset when F-300 key  
is turned

13. **[Main Entrance Door]**

Search ESLA. Open the laser shutter. Assure all laser power supplies  
are ON.

Open the ESLA door.

An audible warning sounds for 5 seconds

The interlock fault LED comes ON

Laser Interlock drops out

The shutter closes

The 3 power supplies remain ON in the MLH

The Excimer power supply in ESLA Turns  
OFF

The YAG power supply in ESLA Turns OFF

Interlock fault will only reset when F-300 key  
is turned

14. **[Roll-Up Door]** Search the Laser ESLA. Open the laser shutter.

Raise the roll-up door up off of the interlock contact.

An audible warning sounds for 5 seconds

The interlock fault LED comes ON

Laser Interlock drops out

The shutter closes

The 3 power supplies remain ON in the MLH

The Excimer power supply in ESLA Turns  
OFF

The YAG power supply in ESLA Turns OFF

Interlock fault will only reset when F-300 key  
is turned

<b>Subject:</b>	<b>SDL End Station Laser Area Interlock Test</b>		
<b>Number:</b>	LS-PPS-0040	<b>Revision:</b>	A
		<b>Effective:</b>	4/13/04
			<b>Page 5 of 6</b>

15. **[Main Door – Door Stop]**

Search the Laser ESLA. Open the laser shutter.

Raise the roll-up door approximately ½” to 1” off of the interlock switch.

An audible warning sounds for 5 seconds

The interlock fault LED comes ON

Laser Interlock drops out

The shutter closes

The 3 power supplies remain ON in the MLH

The Excimer power supply in ESLA Turns OFF

The YAG power supply in ESLA Turns OFF

Interlock fault will only reset when F-300 key is turned

16. **[Reach-back]** Search the Laser ESLA. Open the laser shutter.

Insert a block in the opening of the laser shutter assembly to prevent the shutter from closing.

Open the main entrance door

An audible warning sounds for 5 seconds

The interlock fault LED comes ON

Laser Interlock drops out

The shutter closes

The 3 power supplies turn OFF in the MLH

The Excimer power supply in ESLA Turns OFF

The YAG power supply in ESLA Turns OFF

Reset Fault

17. **[Lockout Switch]**

Secure ESLA for laser & turn On power supply.

Open laser shutter.

Switch from “operate” mode to “Lockout” mode.

Interlock drops out.

The shutter closes

The 3 power supplies remain ON in the MLH

The Excimer power supply in ESLA Turns OFF

The YAG power supply in ESLA Turns OFF

Interlock fault LED comes ON

Attempt to open shutter

Shutter will not Open

Reset interlock fault

Attempt to secure ESLA

ESLA security will not activate

Switch to “Operate” mode

<b>Subject:</b>	<b>SDL End Station Laser Area Interlock Test</b>		
<b>Number:</b>	LS-PPS-0040	<b>Revision:</b>	A
		<b>Effective:</b>	4/13/04
			<b>Page 6 of 6</b>

18. **[Emergency Stops]**

Secure ESLA & turn ON all laser power supplies.

Open laser shutter. Press the external Emergency Stop button (**ES#1**).

The 3 power supplies turn OFF in the MLH

The Excimer power supply in ESLA Turns OFF

The YAG power supply in ESLA Turns OFF

The shutter closes

Interlock drops out in ESLA

Interlock fault is ON

Intermittent audible alarm sounds for 5 seconds.

Interlock fault will only reset when F-300 key is turned

Attempt to secure ESLA.

Laser security will not activate.

Reset emergency stop.

**ES#2**

**ES#3**

19. Secure MLH and ESLA with one person remaining inside & turn ON all laser power supplies.

Press the internal Emergency Stop button.

The 3 power supplies turn OFF in the MLH

The Excimer power supply in ESLA Turns OFF

The YAG power supply in ESLA Turns OFF

Interlock drops out in MLH

Interlock drops out in ESLH

Interlock fault is ON

Interlock fault will only reset when F-300 key is turned

Reset emergency stop.

\* \* \*