

X10A SAFETY CHECKLIST			BROOKHAVEN NATIONAL LABORATORY NATIONAL SYNCHROTRON LIGHT SOURCE
Rev: A	Effective: 1/1/05	Page 1 of 1	Number: LS-SCL-0021
Reviewed by: Thomas McDonald		Reviewed by:	Approved by: A. Ackerman

Original signatures on file.

The only official copy of this file is the one on-line in the NSLS Quality Assurance website. Before using a printed copy, verify that it is the most current version by checking the document effective date on the NSLS QA website.

X10A CAN ONLY OPERATE IN THE FOCUSED MONOCHROMATIC BEAM CONFIGURATION

- ☐ 1. Valid Padlock Index and log, no padlocks open
- ☐ 2. Bremsstrahlung shield (BS #1) in place and banded as per photo
- ☐ 3. Exclusion Zone (EZ #1) wire mesh in place.
 - ☐ A. Beamline Isolation Valve Open
- ☐ 4. Exclusion Zone (EZ #2) wire mesh in place
- ☐ 5. Bremsstrahlung shield (BS #2) in place and banded as per photo
- ☐ 6. Exclusion Zone (EZ #3) wire mesh in place.
- ☐ 7. Bremsstrahlung shield (BS #3) in place and banded as per photo
- ☐ 8. Exclusion Zone (EZ #4) wire mesh in place.
- ☐ 9. Slit view port covered with lead
- ☐ 10. Mirror shields (3) in place, scatter shielding as per photo
- ☐ 11. Exclusion Zone (EZ #5) wire mesh in place
 - ☐ B. Mirror Isolation valve open
- ☐ 12. Bremsstrahlung shield (BS #4) in place and banded as per photo
- ☐ 13. Exclusion Zone (EZ #6) wire mesh in place
- ☐ 14. Exclusion Zone (EZ #7) wire mesh in place
- ☐ 15. Bremsstrahlung shield (BS #5) in place and banded as per photo
- ☐ 16. Exclusion Zone (EZ #8) wire mesh in place
- ☐ 17. Exclusion Zone (EZ #9) wire mesh in place
- ☐ 18. Slit view port covered with lead
 - ☐ C. Monochromator isolation valve open
- ☐ 19. Monochromator leaded glass windows (3) shielded as per photo
- ☐ 20. Exclusion Zone (EZ #10) wire mesh in place
- ☐ 21. Bremsstrahlung shield (BS #6) in place and banded as per photo

Checked by: _____ Date: _____

OPCO: _____ Date: _____