

# **An Overview of Synchrotron X-ray Absorption Spectroscopy**

**Paul Northrup**

Brookhaven National Laboratory  
Environmental Sciences Department  
Environmental Research & Technology Division


NSLS HBCU  
Workshop  
July 19, 2007

“Nothing is so difficult  
but that it may be found out  
by seeking.”  
-Terence (ca. 150 BC)


# Applications of Synchrotron XAS

- A probe of chemical and structural state:
  - Oxidation state and chemical bonding
  - Local and short-range structures
- Element-specific and non-destructive
- Trace or major components, processes
- Applications to chemical, environmental, material and biological systems

# Xray Absorption Spectroscopy:




# XAS measurement




- Direct: transmission through sample.

$$\text{Absorption} = \ln(I_0/I_t)$$

# XAS measurement




- **Indirect**: X-ray fluorescence produced as electron “hole” is filled.
- Characteristic energy for each element.
- Proportional to absorption.




# Three components of XAS:

- Edge step
- Electron transitions
- Extended oscillations
- Each carries different information.




# Absorption edge step




# Fe K absorption edge

- Standards:
  - Hematite:  $\text{Fe}^{3+}$  oxide
  - Vivianite:  $\text{Fe}^{2+}$  phosphate
- Sediment samples at different depths
- Indicative of redox processes





# Electron transitions



- Promotion of electron to available (unfilled) level of absorbing atom -- or neighbor.
- Peak energy differs from edge energy.
- Sensitive to electronic configuration and bonding.

• Rules:

Allowed:

s-p

p-s

p-d

d-p


d-f

Forbidden:

s-s


p-p

d-d




# Uranium L<sub>3</sub> and M<sub>5</sub> edges

- L<sub>3</sub> absorption edge indicates oxidation state
- M<sub>5</sub> edge dominated by 3d > 5f transition
- Importance: U<sup>6+</sup> highly soluble, U<sup>4+</sup> relatively immobile



Rakovan et al



# S K edge

- 2 edge steps (oxidation states)

- 1s to 3p electron transition:

1: sulfide/thiol ( $\text{R-S-R/R-SH}$ ),

2: thiophene,


3: sulfoxide ( $\text{R-(SO)-R}$ ),

4: sulfite/sulfone

( $\text{R-OSO}_2^-/\text{R-(SO}_2\text{)-R}$ ),


5: sulfonate ( $\text{R-SO}_3^-$ ),

6: sulfate ester ( $\text{R-OSO}_3^-$ )




# Organic S species

- Sulfur in sediments
- Sulfate (bio)reduction



- Sulfur in plant roots
- Physiological response to toxin (Zn)





## **“Soft” X-rays:**

- C, N, O edges, very low energy
- Spectral analysis is used to image distribution of different organic compounds and oxides


# EXAFS

- Extended oscillations due to backscatter of electron from neighboring atoms
- Interference pattern:
- Distance
- What element (size)
- Coordination number
- positional/vibrational disorder





# EXAFS data analysis

- Example: S in ZnS
- Convert to frequency (not linear with energy)
- Fourier transform to resolve multiple frequencies



# EXAFS data analysis

S in ZnS structure




Fourier  
transform



# EXAFS data analysis

- Results are obtained by fitting a model to data
- 1D radial distribution, not 3D structure
- Short range



# Limitations

- General:
  - Complexity of spectrum -- all components superimposed
  - Complexity of sample -- shows “Average” configuration
  - Experimental challenges -- sample prep & preservation, energy ranges, interferences, etc.
- EXAFS:
  - requires single species
  - needs greater energy range & concentration
  - fits a model -- model must be close
  - easily abused by over-interpretation