

# **White Paper and Strategic Straw plan for BioXAS Facility at NSLS-II**

Presenter

**Sandeep Rekhi and Mark Chance**

**Case Center for Synchrotron Biosciences at NSLS  
Case Western Reserve University**


# BioXAS related requirements at NSLS-II

- ❑ EXAFS measurements on life sciences samples – existing program
- ❑ High-throughput metalloproteins screening – existing program
- ❑ XAS on screened metalloproteins samples –proposed expansion of the existing program
- ❑ X-ray microprobe studies and imaging microscopy
  - ❑ Time resolved BioXAS


# X3B BioXAS program....future NSLS-II candidate

- ❑ High User Demand (increased to 50% in summer 07, up to 75% in next 5 years)
- ❑ High in-house research usage (metalloproteomics programs and other multidisciplinary projects)
- ❑ High Productivity (~ 45 publication in 2006-2007)
- ❑ High quality service to the users – both technical and scientific
- ❑ Scientific breakthroughs
- ❑ Protein Structure Initiatives (PSI) collaboration – metalloproteomics
- ❑ “Mail-in” and “Drive-by” program
- ❑ Consistency and compatibility between existing X3B at NSLS-I and future BioXAS at NSLS-II
- ❑ State-of-the-art instruments: proposed new monochromator and a 31 element Ge detector (to be moved from X3B)


# Proposed BioXAS beamlines at NSLS-2



# X9B/X3B demand and usage



# X3B/X9B publications



# Significant Upgrades (last five years)

- ❑ New sagittal focusing crystals (achieved ~ 2.5 fold increase in flux)
- ❑ Dynamic focusing using motorized bender (achieved constant spot size)
- ❑ 13 element detector
- ❑ Sample-reference laser alignment
- ❑ New experimental precision bench
- ❑ A state-of-the-art set-up for metalloproteomics program (screen 176 proteins in 7-8 hours)
- ❑ “Mail-in” and “Drive-by” program

# Proposed future upgrades -next 5 years

## ☐ New Monochromator

- Present mono (water cooled) is suitable for a 3-pole wiggler
- New proposed mono will have capability of working at canted DW with both cryo-cooled crystals

## ☐ New Detection system (31 element detector with digital electronics)

- The present detection system can be dedicated to metalloprotein screening (high throughput metalloproteomics program)
- Advantages of new system includes:
  - Low detector saturation (ability to use high flux beam)
  - Ability to measure low concentration samples (increased total active area)
  - Increased data quality (faster data acquisition)
  - High throughput studies (more region of interest selection)


## ☐ Upgrade the metalloproteomics sample holder, stage and sample loading system




# Anticipate additional upgrades by 2015 for smooth transition to NSLS-II

- ☐ New detector and control system
- ☐ Possibility for broader energy coverage
- ☐ Rapid (sub-second) tunability across an edge
- ☐ Planned upgrade of monochromator
- ☐ Cryo sample changer for multiple samples
- ☐ Semi-automatic sample centering system
- ☐ On line sample annealing assembly

# PROPOSED TRANSITION PLAN



# Funding and time Line



# BioXAS at NSLS-II

- ❑ The users time will increase to 80% with 20% for developmental activity
- ❑ Expansion of in-house metalloproteomics program
- ❑ Full support from the CSB team to run the life science XAS program, help users to submit proposals and beamtime management
- ❑ Expansion of mail-in and drive-by program

END

**Slide 1:** Summary of long term scientific program at X9B/X3B, Impact on science, Publications, Significant upgrades, Proposed upgrades in next five years,

**Phase 1** Proposed transfer of X3B hutch components to X3A and X3B beamline to a 3-pole wiggler line

**Phase 2a** Proposed transfer of X3A beamline to a 3-pole wiggler line

**Phase 2b** Cooperation with the DW project specialized beamlines for the cutting edge life science related XAS experiments such as flow cell, time resolved XAS experiments, microprobe XAS experiments

**Phase 3** Building an additional canted wiggler XAS beamline in conjunction to the proposed DW footprinting beamline

**Slide 2:** Transition plan:

Option 1: Move X3B to X3A while transition

**Slide 3:** Staff support in running the life science program, helping users to submit proposals and beamtime management, etc.

Slide 4: Global impact on users → 80% time for user community and 20% for developmental activity


Growth of in-house projects such as metalloproteomics

Transfer will require smoother and quicker transition

# Metalloproteomics in Structure Genomics

- ❑ High-throughput characterization of metalloproteins using X-ray absorption spectroscopy
- ❑ Developed in structure genomics to assist crystallographic phasing, functional and structural annotation of metalloproteins
- ❑ 2164 (654 in PSI – 1 and 1510 in PSI – 2 NYSGXRC proteins analyzed and 10% showed presence of a transition metal (Zn, Mn, Ni, Co and Cu)
- ❑ 95% accurate based on 143 solved crystal structures
- ❑ Long term goal: to establish a metalloprotein database  
~ 1000 families are covered so far

## Experimental Setup at Beamline X3B



# Lab Space Specification and Location

- ☐ Clean Room
- ☐ Fume Hood
- ☐ Oven and Furnace
- ☐ 80 K deep freeze refrigerator
- ☐ Centrifuge and related instruments
- ☐ Sample storage and cold room
- ☐ Sample preparation work bench space
- ☐ Storage space for excess equipments
- ☐ Electronics and beamline instrumentation construction lab

## **Office Space Specification and Location:**

- ☐ Near to the beamline if possible
- ☐ 1 or 2 staffs per office room

## **Instrumentation Lab Space and Location :**

- ☐ Near to the beamline
- ☐ Large enough to provide a work space of at least 3 personnel



## **Key attractions of XAS for Biological Research:**

- (1) XAS provides information on electronic and atomic structure for both crystalline and non-crystalline systems;
- (2) XAS is a comprehensive atomic level structural tool sensitive to within  $\sim 4\text{-}5\text{ \AA}$  of metal sites;
- (3) XAS allows an order of magnitude more accurate bond length determination than that obtained by protein crystallography;
- (4) XAS is an extremely fast probe ( $\tau < 10\text{-}14\text{ s}$  for Fe *K-edge*) *that is* suitable for multiple-scale time-resolved experiments; Combination of (1)&(4) makes XAS a unique technique for probing structure of reactive intermediates in solutions.
- (5) Contrary to UV-VIS and EPR spectroscopy XAS is always detectable: There are no “spectroscopically quiet” metals;
- (6) XAS is capable of probing dilute samples at the micromolar level;
- (7) Since XAS measurements are usually done at 10-20K, biosamples are less susceptible to x-ray beam photoreduction compared to the XRD. Photoreduction at the active center is easy to monitor via XANES.