



actinoids

- New materials and smaller features drive need for characterization techniques
- Need compositional and electrical information at device dimensions

IBM at NSLS-I

■ *in-situ* real-time annealing effects

- ❑ phase formation in silicides
- ❑ diffusion barrier failure T
- ❑ solder formation
- ❑ crystallization of phase-change materials




■ phase identification (GIXD)

- ❑ high-κ dielectric ultra-thin films, intermetallics, phase-change materials, silicides

■ texture and grain size (pole figures, GIXD, HRXRD)

- ❑ axiotaxy in silicides and germanides, Cu features (interconnects)




■ surface and interface quality (XRR)


- ❑ dielectrics, silicides, phase-change materials, etc.

■ strain in thin films and features (HRXRD, microdiffraction)


- ❑ SiGe/Si, Si:C/Si, etc.



*to understand and
improve materials and
processes used in
making devices*




Endstation instrumentation at X20




"Standard" diffraction

- Standard diffraction
- Triple-axis diffraction
- Reflectivity
- GIXD
- SAXS
- Four-circle Huber used for pole figures, reflectivity, structure ID, strain relaxation, etc.



Time-resolved scattering



- High temp. annealing chamber
- Diffraction into fast linear or area detector
- Resistance and optical scattering
- Used with multilayer monochromator for *in-situ* measurements of phase transformations, barrier failures, kinetics, texture evolution.



Microdiffraction

- Standard diffraction
- Fluorescence
- Heating/electrical stage
- Used for diffraction topography, grain or local tilt and mosaic, elemental mapping.

Time-resolved *in-situ* XRD, resistance and light scattering




- **Fast data collection and analysis**
- **Intense signal**
 - Transformation temps
 - Morphology changes
 - Kinetics
- **Quick evaluation of large array of variables**
 - processes parameters
 - film thicknesses
 - dopants
- **Size effects**

limitations:

- **detector**
- **reciprocal space access**
- **sample turnaround**
- **?**


What will microelectronics look like in 8 years?





• Gate Oxide no longer scales:

- Need to reduce gate leakage
- Power Density increases – more transistors with constant power
- Power constant – Gate can not turn off channel effectively if gate oxide is increased

Fundamental Atomic Limit to Scaling Recipe



Oxide thickness is approaching a few atomic layers



?

As long as IBM makes hardware, we'll need materials and process analysis...

IBM Needs at NSLS-II

- High-throughput/rapid access very important:
 - **phase ID, texture, failures, phase formation, process variables**
- **Guaranteed** access to high-throughput time/temp-resolved XRD
- Access to thin film analysis techniques
 - GIXD, texture, XRR, Bragg-Brentano, etc.
- Access to micro/nano diffraction
- Possible access to:
 - nanotomography (Xradia full-field microscope?)
 - EXAFS
 - PDF?
- Minimum time gap in access between NSLS and NSLS-II

Tentative plans

- Upgrade of **high-throughput** TR endstation on X20C (in stages):
 - ❑ dome for more reciprocal space access
 - ❑ special diffractometer
 - ❑ RTA heater and laser annealing capabilities
 - ❑ ports for adding probes
 - ❑ area, linear and point detectors
 - ❑ automated sample changing
 - ❑ automated evacuation and ambient control
 - ❑ remote access
- New 6-circle for X20A
- Move both to NSLS-II
- If you are interested, send email to JLJ@BNL.GOV