

Parameters affecting crystal lifetime in MX and possible radiation damage mitigation strategies.



**MX Frontiers at the
One Micron Scale
BNL, USA.
July 2009**






Elspeth Garman, LMB, Oxford
Elspeth.garman@bioch.ox.ac.uk
Karthik Paithanker, LMB, Oxford

The Plan:

- **A metric for Radiation Damage.
Dose: RADDOSE.**
- Scavengers: RT and 100K.
- Simultaneous multi-crystal data collection and data retrieval.

What really happens when X-ray photons hit the crystal ?



$\lambda = 1 \text{ \AA}$ (at energy 12.4 keV) for a
100x100x100 μm crystal

Compton scattering
and photoelectric effect both
deposit energy in the crystal

DOSE

- DOSE is the ENERGY lost per KILOGRAMME (!!)
- Measured in Joules/kg i.e. the absorbed energy per unit mass.
- Fundamental metric against which to measure damage.
- FLUX is in photons/second.
- Flux density is in photons/second/unit area.
- Takes care of the physics but NOT the chemistry.

DOSE Postulate:

- Damage at 100K is proportional to dose (I_{mean} , Bfactor, R_{merge} , R_d , specific structural damage [dose rate effect?]).
- There is a MAXIMUM dose (Joules/kg = Gy) which protein crystals can tolerate which depends only on the PHYSICS of the situation.
- Crystal might not reach that limit due to chemical factors, but it will not last BEYOND the limit.
- Need to be able to calculate this DOSE:
[RADDPOSE: Murray, Garman & Ravelli, JAPC 2004]

Way of estimating absorbed dose,

$$\mathbf{D: (Gy = J\ kg^{-1})}$$

Dose rate = mass absorption coeff * photon energy *
number of photons in unit time / Area

$$\mathbf{dD/dt = (\mu/\rho)\ E\ I_{inc}} \quad (\mathbf{I_{inc}} = \text{incident flux density})$$

For (μ/ρ) in cm^2/g , I_{inc} in $\text{photons/s}/\mu\text{m}^2$, E in keV, t in seconds,
total dose is:

$$\mathbf{D = (\mu/\rho)\ E\ I_{inc}\ t\ 10^{11}\ (Gy)}$$

e.g. $(\mu/\rho) = 2.6\ \text{cm}^2/\text{g}$ (50% solvent), $E = 12\ \text{keV}$,
 $t = 100\ \mu\text{s}$, $A = 80 \times 80\ \mu\text{m}^2$ beam cross section


$$\mathbf{D = 7.8 \times 10^{-8}\ Gy/photon}$$

- For $10^6\ \text{Gy}$, 1 ionisation / 20 amino acids for a 400 a.a. protein molecule. [See O'Neill, Stevens & Garman. JSR (2002) **9**, 329-332]

Can define Decay R_{factor}
to plot against $\Delta(\text{dose})$:

R_d : pair wise R factor between identical and symmetry
related reflections occurring on different diffraction
images.


$$R_d = \frac{\sum_{hkl} \sum_{|i-j|=d} |Y_i - Y_j|}{\sum_{hkl} \sum_{|i-j|=d} (Y_i + Y_j)/2}$$



Can define a s_{AD} (to plot against dose):

- Coefficient of sensitivity \propto change in relative isotropic B factor: [Kmetko et al 2006, Acta D62, 1030]

$$s_{AD} = \Delta B_{rel} / 8\pi^2 \Delta D \quad (\text{e.g. HEWL} = 0.012 \text{ at } 100\text{K})$$



Make this easier for MX
(include solvent contribution in mM
and heavy atoms explicitly)


RADDOSE

Use crystal and beam characteristics to calculate
the dose.

Calculating Dose


(*RADDOSE*)

Beam Characteristics



Beam Characteristics


Beam profile




ID14-4

Wavelength




Photons per second



Calibration



Beam Characteristics






Pin-diode current vs.
integrated counts at different
attenuations
Hamamatsu Si (S3204-09)

Photons per second vs.
integrated counts

N.B.; Two calibrations required.

Determination of X-ray flux using silicon pin diodes. RL Owen, JM Holton, C Schulze-Briesse, EF Garman. *JSR* (2009) **16**, 143-151.

Beam Profile Comparison uniform or 2D Gaussians




Beam Profiles may be entered into RADDOSE as uniform or 2D Gaussians

- Lifetime curve is reciprocal of intensity curve
- Differential irradiation may lead to differential damage

get data which merge poorly

- Crystal heating might be higher than first predicted




Number of amino acids

‘HA’ atoms per monomer, *e.g.* S, Se, Hg

Solvent - concentrations of components, *e.g.* Na⁺, Cl⁻

Calculating Dose (*RADDOSE*)

Crystal Characteristics




absorption coefficients at 12.4 keV

e.g. apoferritin: 0.406mm^{-1}

holoferritin: 1.133mm^{-1}

Experimental Absorption Coefficients for heavy atoms

Variety of anomalous edges




f'' is proportional to μ_{pe}

- One can normalise a fluorescence spectrum to known values of
- μ_{pe} far from the absorption edge.
- This is implemented in RADDOSE using the SPLINOR file from CHOOCH. (G. Evans and R. F. Pettifer *J. Appl. Cryst.* **34**, 82-86, 2001.)
- Other absorption values are taken from library values McMaster 1960 and mucal.f

Calculating Dose (*RADDOSE*)

Crystal Characteristics

Beam Characteristics



Experimental Dose Limit (100K)

For $I_0 = 1/2$

$$D_{1/2} = 4.3 (\pm 0.4) \times 10^7 \text{ Gy} = \mathbf{43 \text{ MGy}}$$

(cf 'Henderson limit' $20 \text{ MGy} \equiv 5 \text{ electrons}/\text{\AA}^2$)

$43 \text{ MGy} \cong 10 \text{ electrons}/\text{\AA}^2$

cf. hamster death 3 Gy)

Suggested limit to retain biological 'fidelity'

$I_0 = 0.7$

$$D_{0.7} = 3.0 \times 10^7 \text{ Gy} = \mathbf{30 \text{ MGy}}$$

$D_{0.7}$ for ferritin corresponds 10^7 photons/unit cell

Robin Leslie Owen, Enrique Rudiño-Piñera, Elspeth F. Garman.

PNAS (2006) 103, 4912 - 4917.

Assumptions in distributed version of RADDOSE

- $\mu_{\text{abs}} = \mu_{\text{pe}}$ small underestimate at high energies: as Compton scattering is neglected
- Fluorescent X-rays are absorbed:
results in an overestimate of dose for heavy scatterers
- Crystal rotation neglected
- No potential dose rate effect considered.

‘X-ray Absorption by Macromolecular Crystals; the Effects of Wavelength and Crystal Composition on Absorbed Dose’. Murray, Garman, Ravelli, *J. Appl. Cryst.* (2004) 26, 513-522

Improved RADDOSE (2009)


- Input/Output made independent of CCP4 libraries
- Outputs time to reach experimental dose limit (30 MGy)
- Makes a correction related to physics of energy loss

For atoms with $Z > 20$ after interaction via photo electric effect the electron can relax via

1. Auger effect
2. X-ray fluorescence



This fluorescence energy may escape from crystal.

- Energy loss from Compton scattering included.





Why care about X-ray fluorescence escape ?

- If the incident energy is **greater** than the absorption edge energy; that atom may undergo photoelectric excitation
- Atom can decay *via* Auger or X-ray fluorescence
- X-ray fluorescence photon can escape, depending on its energy, the thickness of crystal, decreasing the energy lost in the sample.
- **May be important for micro crystals**
- Knowledge of correct dose \rightarrow correct estimation of lifetime \rightarrow planning of the experiment



X-ray fluorescence escape:

- For C, N, O, S there is low probability of X-ray fluorescence.
- For heavy elements ($Z > 20$) such as Se there is higher chance of X-ray fluorescence
- The contribution to the anomalous scattering factor (f'') is directly proportional to the photoelectric cross-section.
- **Corrected energy**
= deposited energy
 - K -shell escape, – L_I escape,
 - L_{II} escape, – L_{III} escape




e.g. SeMet protein – phosphomethylpyrimidine kinase

- Seleno-methionine crystals often have a shorter lifetime than native crystal, due to the large photoelectric cross section of selenium.
- Using RADDOSE can predict the maximum crystal lifetime for a MAD experiment: predicted lifetime increases by 27 % at 12.6634 keV when the fluorescent escape is included.


Crystal size:
0.04 0.1 0.05 mm³

Tophat beam:
0.1 0.1 mm²




Inclusion of energy absorption due to Compton scattering

- Compton scattering: interaction of a photon with a 'free' electron (outer shell)
- The incoming photon will inelastically scatter and the electron will recoil at an angle φ
- The energy of the recoil electron is deposited in the crystal itself (contributes to the absorbed dose)
- Total energy loss $\propto \sigma_{\text{Compton cross-section}} E_{\text{recoil electron}} + \sigma_{\text{photoelectric cross-section}} E_{\text{photoelectric absorption}}$




Maximum electron recoil energy:

$$E_{\text{recoiled electron}} = \frac{2E_{\text{incident}}^2}{mc^2 \left(1 + 2\frac{E_{\text{incident}}}{mc^2}\right)}$$



Effect of Compton scattering energy loss on dose



PPK, 0.04 0.1 0.05 mm³ crystal, tophat beam, 0.05 0.05 mm²,
10¹² photons/mm², 0.2 s/image. 398 residues. 12 seleniums

Diffraction-dose efficiency

- Want to maximise:

$$\frac{\text{Diffracted intensity}}{\text{Dose}} = I_{\text{DE}}$$


$$\frac{I_{\text{scatt}}}{\text{Dose}} = \frac{I_{\text{scatt}} \times V \times \rho}{\text{Energy absorbed}} \propto \frac{\lambda^2 e^{-\mu t} \times V^2 \times \rho}{\text{Energy absorbed}} \quad \text{N.B.}$$

[Dose=Energy absorbed/mass]

$$\frac{I_{\text{scatt}}}{\text{dose}} \propto \frac{\text{Volume}^2 \times \lambda^2 \times e^{-\mu_{\text{att}} \text{thickness}}}{en_{\text{incident}} (1 - e^{-\mu_{\text{photoelectric}} \text{thickness}}) + en_{\text{Comptonelectron}} (1 - e^{-\mu_{\text{Compton}} \text{thickness}})}$$

I_{DE} : Is there an optimal energy for a given sized sample ?


Protein:
Phosphomethyl
pyrimidine
kinase P1
(398 residues)
5 sulphur atoms
12 seleniums
Beam: tophat
 10^{12} photons/sec
 $x, y = 0.3 \times 0.3$
 mm^2



I_{DE} : effect of selenium for a 100 μ m thick crystal

100 micron x 100 micron beam; 1E12 photons/sec; 398 residues


100 micron x 100 micron x 100 micron crystal



IDE: effect of selenium for a 1 μm crystal


1 micron x 1 micron beam; 1E12 photons/sec; 398 residues

1 micron x 1 micron x 1 micron crystal



I_{DE} : effect of heavy atoms: 500 μm crystal

376 residues, 8 monomers in unit cell, 3 Uranium atoms per monomer in derivative
47% solvent content, beam=crystal size 0.5 x 0.5 x 0.5 mm³




human phosphatase binding protein (HPBP)

I_{DE} for a 5 μm crystal.

376 residues, 8 monomers in unit cell, 3 Uranium atoms per monomer in derivative


47% solvent content, beam=crystal size 0.005 x 0.005 x 0.005 mm³



human phosphatase binding protein (HPBP)


Crystal and beam size corrections

- ◆ Assume blue box is the X-ray beam perpendicular to the screen and black arrow is the rotation axis
- ◆ For crystals whose dimensions exceed that of the beam, dose as calculated for a stationary crystal is not an accurate metric for the estimation of radiation damage
- ◆ Taking into account the irradiated volume of the crystal if the crystal is bigger than the beam
- ◆ Knowledge of physical orientation of the crystal with respect to X-ray beam




The Plan:

- A metric for Radiation Damage. Dose: RADDDOSE.
- **Scavengers: RT and 100K.**
- Simultaneous multi-crystal data collection and data retrieval.



The chemistry:

mobile e^- go to electron affinic sites





Specific structural damage

DISULPHIDE BONDS (S-S) MOST SUSCEPTIBLE

Weik *et al* (2000) PNAS 97, 623-628

Burmeister (2000), Acta Cryst D56, 328-341.


Ravelli and McSweeney, (2000) Structure 8, 315-328



Water/solvent chemistry






1) Ionization

Ionizing radiation



2) Electronic excitation




Ionizing radiation






Scavengers at cryotemperatures:

Rationale


- Crystals are usually cooled to 100K to reduce the mobility of free radicals.
- The existence of specific damage at this temperature shows that some species are still mobile (electrons).
- Therefore scavengers may be able to react with these species and reduce their mobility and the reactivity of fixed species, protecting the crystal from specific damage.






Cystine




Cystine



Ascorbate




Cystine




Quinone


(a) HEWL
native dataset 1




(b) HEWL
native dataset 6



(c) HEWL
ascorbate
dataset 1



(d) HEWL
ascorbate
dataset 6



Change in atomic B factors of refined structures with dose.

No increase in temperature factor




70 % increase in temperature factor

(a) HEWL
native dataset 1



(a) HEWL
native dataset 6




(a) HEWL
ascorbate dataset 1




(a) HEWL
ascorbate dataset 6




Potential radioprotectants identified by on-line microspectrophotometry



Ascorbate




1,4-Benzoquinone



2,2,6,6-tetramethyl-4-piperidone (TEMP).

RT: Ascorbate, co-crystallised 1M



$$D_{1/2} \text{ Native} = 0.125 \text{ MGy}$$

$$D_{1/2} \text{ Ascorbate} = 0.24 \text{ MGy}$$

Native–Exponential (1st Order), Ascorbate–Linear (0th Order).

Decay of native crystals is linear at 100K – is the RT native exponential decay dominated by OH radicals?

Dose rates (Gy/Sec) = Native 1 – 6, Native 2 – 12.8

Ascorbate 1 – 6.4, Ascorbate 2 – 6.4


RT: ascorbate ESRF data: 2800 Gy/s

Synchrotron Data


$$D_{1/2} \text{ Native} = 0.9 \text{ MGy (x8)}$$

$$D_{1/2} \text{ Ascorbate} = 2.2 \text{ MGy (x9)}$$

$$D_e = 2.$$



Quinone, soak, 14 days, 1M



$D_{1/2}$ Native = 0.14 MGy,

$D_{1/2}$ Quinone = 1.2 MGy

$D_e = 9.$

Native – Exponential (1st order), Quinone – Linear (0th).

Dose rates (Gy/Sec) = Native 1 – 6.4, Native 2 – 6.0

Quinone 1 – 5.7, Quinone 2 – 5.9

Electron density difference map analysis shows no specific damage.

Scavengers

- Effective at RT (benzo-quinone factor of 9, ascorbate factor of 2)
- Higher dose tolerance when using scavengers with higher k_{e-} .
(Barker, A.I., Southworth-Davies, R. J., Paithankar, K.S., Carmichael, I. and Garman, E.F. *J. Synchrotron Rad.* (2009). **16**, 205–216)
- $De = D1/2(\text{scavenger}) / D1/2(\text{native})$


Sodium Ascorbate $k_{e-} = 3.0 \times 10^8 \text{ M}^{-1}\text{s}^{-1}$ $k_{OH} = 8 \times 10^9 \text{ M}^{-1}\text{s}^{-1}$ $De = 2$

1,4-Benzoquinone $k_{e-} = 1.2 \times 10^{10} \text{ M}^{-1}\text{s}^{-1}$ $k_{OH} = 1.2 \times 10^9 \text{ M}^{-1}\text{s}^{-1}$ $De = 8.9$

- Not clear yet at 100K. Ascorbate is effective in suppressing specific damage.
- Is it possible to isolate $e^-(aq)$ and OH^\bullet effects?
- Try an electron scavenger:


Sodium Nitrate $k_{e-} = 9.7 \times 10^9 \text{ M}^{-1}\text{s}^{-1}$

- 100K: Lysozyme solution (NaAc 200mM pH 4.7, 10%w/v NaCl lysozyme 50 mg/mL, 20 % glycerol)
- $e^{-}(aq)$ produced during the first fs of irradiation produced the disulfide radical anion $Cys\cdot$ with a characteristic peak around 400 nm.



100K: Lysozyme solution + 1M
sodium nitrate

$e^-_{(aq)}$ scavenged effectively.
No S-S \cdot are observed.




R_d

Lysozyme crystal + 1M sodium nitrate

$D_e = 3$

Protects disulfide bonds up to 70 MGy.



Radioprotectants: Conclusions

Not yet seen more than a factor of 3 in global damage at 100K, but have seen protection of amino acids so definitely worth considering for specific cases.

Have potential to make a significant difference at room temperature.

We are working on trying to understand the changes from first order to zeroth order kinetics at RT.




The Plan:

- A metric for Radiation Damage. Dose: RADDOSE.
- Scavengers: RT and 100K.
- **Simultaneous multi-crystal data collection and data retrieval.**

Objectives

Total Cryst

- Adapt methods from material sciences to macromolecular crystallography
- Use of multiple crystals is routine in chemical crystallography
- Take multiple crystals (2, 3, ... n) in a single loop and collect data – oriented randomly
- Index the diffraction pattern and utilize the information from all the crystals
- **Test the feasibility** of the new methods for MX to combat radiation damage
- Computationally provide accurate estimate of the maximum dose limit using the program RADDOSE
- Optimise the energy incident for a given crystal size, composition



Why ? Specific structural damage

Radiation damage in crystallography even at 100 K



Usual: Collect single crystal dataset for 100

Proposed: collect 20 of data from 5 different crystals simultaneously and combine them

data of 20 5 crystals ~ data of 100




Advantages:

- ✓ All the 5 crystals are at the beginning of decay
- ✓ Lower absorbed dose per crystal
- ✓ Higher quality data
- ✓ Metal centres, active site preserved
- ✓ Extract correct biological information




Potential applications of multiple crystal diffraction

- ◆ Micro crystals used in structure solution
- ◆ Minimize crystal handling (less mechanical damage)
- ◆ Multiple micro crystals could be mounted with Crystal catcher system
- ◆ Animal hairs attached or woven (like triple helices) on the surface of loops could be used to fish micro-crystals
- ◆ Streak seeding could be done with such loops and subsequently left in the drop (protein + precipitant) to grow crystals on their surface




Crystal Catcher: Kitatani *et al.*, (2008) Appl. Phys. Express **1**, 370021-3

Multi-crystals: abundance of reflections leads to spot overlap





Overlap fraction with multiple crystals (mosaicity = 0.5°)


Dotted lines → Lysozyme
Solid lines → Insulin



Simulations of diffraction patterns




Diffraction patterns with 1, 2, 3 and 4 crystals of lysozyme (top) and insulin (bottom), $\Delta\phi = 0.5^\circ$






Experiments: Multiple crystals in X-ray beam


Data of crystals from
lysozyme and insulin
collected at ESRF
(ID 14-4)



Insulin




Lysozyme




Data reduction in practice using materials science software

- Search for peaks above a certain threshold and construct an array of peaks
- Indexing with pattern recognition algorithm - Grainspotter
- **OUTPUT:** orientation of each crystal in the ensemble
- The orientation of each crystal lattice is provided to **MOSFLM** by means of the U (orientation) matrix and data integrated
- Red and green indicate the two different lattices from two crystal lysozyme dataset




www.totalcryst.dk (all programs released as open source: www.sf.net)

A single lattice from a four crystal Lysozyme dataset



Animation
illustrates the
four lattices
identified
uniquely from
a dataset
obtained from
a multiple
crystal data
collection
with four
crystals in the
beam



Summary – *TotalCryst*

Simulations show spot overlap is not huge problem; due to random overlap and not systematic overlap

Data were collected with multiple single crystals in a single loop

Possible to index unknown lattice but collecting few exposures with only one crystal or select an exposure with one strong lattice

Extraction of data from up to 7 crystals achieved

The combination of data from multiple crystals compensates for the loss of redundancy owing to rejected spots

Need to extend the experiments to microcrystals.

Paithankar, K.S., Sørensen, H.Ø., Poulsen, H. F & Garman, E. F. *in preparation*

Thanks to:

My group, past and present:

Karthik Paithankar
Eugenio de la Mora Lugo
Adam Barker
James Murray (IC)
Robin Owen (DLS)
Enrique Rudiño-Piñera (UNAM)
Robert Southworth-Davies (DLS)



Raimond Ravelli



Martin Weik

Ian Carmichael (Notre Dame)
John McGeehan (U Port)
Sean McSweeney (ESRF)
Raimond Ravelli (ULMC)
Martin Weik (IBS)
Henning Sorensen (Riso)
Jon Wright (ESRF)



**Adam
Barker**



Ian Carmichael

**ESRF Grenoble for LS2047,
MX-161, MX348, MX438, MX-666,
MX-812, MX-931**

Funding:



Total Cryst
EU FP6



Oxford Cryosystems

