X1A

J.A. Brandes, G.D. Cody (CIW), J.I. Hedges (U. of Washington), and S.Wirick (SUNY at Stony Brook)

The structural and molecular composition of particulate organic carbon was examined using the imaging and spectral capabilities of the X1-A STXM. The oceanic carbon cycle starts with photosynthetically produced organic matter in surface waters, which, after death of the organisms, settles to the ocean floor. During this settling process, roughly 50-90% of the organic matter is remineralized to CO₂. Although much work has been done to characterize the bulk chemical changes in this particulate flux as it settles, the physical structural changes have been difficult to elucidate. Soft X-ray microscopy and C-XANES spectroscopy has applied to a series of samples collected in the Equatorial Pacific ocean at depths between 105 meters and 3450 meters. Figure 1 shows a sample from 3450 m depth that exhibits regions of both highly degraded and pristine cellular material, including identifiable cell wall structures. The different carbon types are located within a few microns of each other, and a region can be observed in the lower right side of Figure 1 where cellular wall structures convert within a distance of less than a micron to more degraded carbon. These features are located together with a matrix of calcium carbonate and silicate phytoliths. Images taken at different energies permit the identification and location of these different components within the complex sample matrix.




Figure 1.