

Brookhaven National Laboratory National Synchrotron Light Source		Number: LS-PPS-0027	Revision: A
		Effective: 6/14/02	Page 1 of 3
Subject: <u>X14B</u> Radiological Interlock Test			
Prepared/ M. Buckley		Approved By: S.Buda	
Approved By:			

*Approval signatures on file with master copy.

[Revision/Periodic Review Log](#)

Test Reason:	Test Result: <input type="checkbox"/> Passed <input type="checkbox"/> Failed	
	Test Type: <input type="checkbox"/> Full <input type="checkbox"/> Partial	
Test Date:	Start Time:	Finish Time:
Tester 1:	Assistant 1:	
Tester 2:	Assistant 2:	

Preparation:

Inform control room operator that test will be done.

Turn off 'Red Tag' switch on x-ray ring interlock, and apply LOTO to X-Ray Security system.

Verify that vacuum valves and water interlocks are OK on the beam lines that will be checked for Phase I beamlines or connect user interlock test jumper.

- Place latch device on door lock so key can be removed. Place switch holders on door switches. Search hutch, put door key in SRU and turn and open safety shutter.

The "Beam On" sign turns ON

Remove holder from switch 2.

RIB latch light comes ON

Hutch	
Top	Bottom

- Replace holder on switch 2.

RIB cannot be reset unless SOR key is turned.

- Remove holder on switch 1.

Safety shutter closes
RIA latch light comes on
RIA can't be reset w/o SOR

- Replace holder on switch 1. Interlock hutch, leave shutter closed.

Hutch Kirk key (KKH) cannot be removed from SRU unless button is pushed

- Open safety shutter.

KKH cannot be removed even if button is pushed

- Close shutters. Turn KKH and note where solenoid unit stops rotation. Push the button and turn key just beyond this point, but not far enough to actuate "key removed" switch. Release the button and leave key in this position. Open shutter, remove KKH, and replace.

RIA latch light comes on
RIB latch light comes on
The safety shutter closes

X14B Radiological Interlock Test			
Number: LS-PPS-0027	Revision: A	Effective: 6/14/02	Page 2 of 3

7. With hutch interlocked, leave shutter closed.

Beam line transfer key can be removed only if
button on SRU is pushed _____

Open safety shutter.

Beam line transfer key can't be removed even if
button is pushed _____

8. Close safety shutter. "Cheat" SRU as described in step 6. Open safety shutter, remove the key and replace.

RIA latch light comes on _____

RIB latch light comes on _____

The safety shutter closes _____

SPA comes on _____

SPA can't be reset w/o SOR _____

Reset RIA and RIB

9. Open the safety shutter. Remove the Control Room Shutter Enable Key for this beam line, then replace the key.

The safety shutter closes _____

RIA latch light comes on _____

The "Control Room Shutter Enable" light goes out _____

The SPA latch comes on (after the key is back on) _____

The safety shutter cannot be opened _____

Reset RIA and SPA

10. Open the Safety Shutter.

Turn "Auto Open" key to "ON".

At the Control Room Operator's Console, switch the shutter master enable to "Inhibit" and back to "Enable"

The safety shutter closes _____

And then re-opens _____

Turn "Auto Open" key to "OFF".

Switch the master enable to "Inhibit" and back to "Enable"

The safety shutter closes _____

And does NOT re-open _____

X14B Radiological Interlock Test			
Number: LS-PPS-0027	Revision: A	Effective: 6/14/02	Page 3 of 3

11. One person goes to the x-ray equipment area and resets the indicators for RIAX and RIBX at SR100.

Open the safety shutter.

Remove switch holder from door switch number 2.

RIB latch light comes on _____

RIBX "Loop Enabled" indicator goes out _____

RIBX "Loop Disabled" indicator comes on _____

Replace switch holder.

RIBX cannot be reset _____

Reset RIB

RIBX CAN be reset

("Loop Enabled" indicator comes on) _____

12. Remove holder from door switch 1.

RIA latch light comes on _____

RIAX "Loop Enabled" indicator goes out _____

RIAX "Loop Disabled" indicator comes on _____

The safety shutter closes _____

RIAX cannot be reset _____

Reset RIA

RIAX CAN be reset

("Loop Enabled" indicator comes on) _____

13. Remove the holders from the switches. Remove the latch device from the door lock(s). Stand inside the hutch and open and close the door.

Both door switches operate freely and each makes a "click" when door is opened

(Complete test for both door switches)

Hutch Door Switch	
Top	Bottom

14. Search the hutch. Locate the safety shutter associated with the beam line by tracing the beam pipe through the shield wall. Do not rely on labels on the shutter mechanism. Open the shutter.

The correct shutter opens _____

No other shutter at that saw tooth opens or attempts to open _____

Close the safety shutter.

The mechanism moves freely and without hesitation _____

15. Remove user interlock test jumper (if applicable). _____

* * *

The only official copy of this file is the one on-line in the NSLS Quality Assurance website. Before using a printed copy, verify that it is the most current version by checking the document effective date on the NSLS QA website.